


Intelligent

?LR@5MRGA



Classification & Composition

Classification method of intelligent instantaneous sterilizing equipment



Aj_qqj_a_rgmI_aampbjle rm
bj_cpcIrfc_rgle k crfmbq



Bjpcar_fc_rgle



Slbjpcar_fc_rgle



Aj_qqj_a_rgmI_aampbjle rm
bj_cpcIrbjgldcargmI k crfmbq

Qrc_k_bjgqldcargmI
kqrc_k_bjgqldcargmI'

Qsncpfc_rcb_u_rcp_bjgqldca+
rgmI ku_rcp_bjgqldcargmI'

Intelligent instantaneous sterilizing equipment is divided into 4 series

- +
- +
- +
- +



Rwnc 6 gIrcjjjecIrgIqr_lrcIcmsq qrcpjpxgle cosjn k cI r+ bjpcar
fc_rgle_l b qrc_k_bjgqldcargmI

Rwnc 66 gIrcjjjecIrgIqr_lrcIcmsq qrcpjpxgle cosjn k cI r+
bjpcar_fc_rgle_l b u_rcp_bjgqldcargmI

Rwnc 666 gIrcjjjecIrgIqr_lrcIcmsq qrcpjpxgle cosjn k cI r+
gIbjpcar_fc_rgle_l b qrc_k_bjgqldcargmI

Rwnc 6T gIrcjjjecIrgIqr_lrcIcmsq qrcpjpxgle cosjn k cI r+
gIbjpcar_fc_rgle_l b u_rcp_bjgqldcargmI

K_gIrf_gIcp



Q5+/0...NJA

/...k q

MNA BN KMB@SQ
QCPTCP AJXCLR

MB@A-QQJ

BAQ

